



Product Datasheet

VAV3 Antibody (orb1246802)

Description	VAV3 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: C-SGEQGTKLPEK	orb1246802 (0.5 ug/ml) staining of Human...
Target	VAV3	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1246802 (7.5 ug/ml) staining of paraf...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 100kDa band observed in Human Bone Marrow lysates (calculated MW of 97.8kDa according to NP_006104.4). Recommended concentration: 0.5-1.5ug/ml. Primary incubation 1 hour at room temperature. Immunohistochemistry: Paraffin embedded Human Spleen. Recommended concentration: 7.5ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q9UKW4	orb1246802 Negative Control showing stai...
NCBI	NP_006104.4 , NP_001073343.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 100kDa band observed in Human Bone Marrow lysates (calculated MW of 97.8kDa according to NP_006104.4). Recommended concentration: 0.5-1.5ug/ml. Primary incubation 1 hour at room temperature. Immunohistochemistry: Paraffin embedded Human Spleen. Recommended concentration: 7.5ug/ml.	
Expiration Date	12 months from date of receipt.	