

Product Datasheet

**IL18 Receptor Beta antibody (PE) [Out of stock]
(orb124683)**

Description

Rabbit polyclonal to IL18 Receptor

Species/Host	Rabbit
Conjugation	PE
Immunogen	KLH conjugated synthetic peptide derived from human IL-18R beta
Target	IL18RAP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Application notes	orb124683 is the PE conjugated form of orb5545
Isotype	IgG
Clonality	Polyclonal
MW	68kDa
Uniprot ID	Q9Z2B1
Expiration Date	12 months from date of receipt.