


Product Datasheet

TIRAP Antibody (orb1246849)

Description		
TIRAP Antibody		
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	orb1246849 (0.1 ug/ml) staining of K562 ...
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-QTLLKKPKKRPNSPE	
Target	TIRAP	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 24kDa band observed in lysates of Human myeloid erythroblast call line K562 (calculated MW of 23.9kDa according to NP_001034750.1 and 24.8kDa according to NP_683708.1). Recommended concentration: 0.1-0.3ug/ml.	
Clonality	Polyclonal	
Uniprot ID	P58753	
NCBI	NP_683708.1, NP_001034750.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 24kDa band observed in lysates of Human myeloid erythroblast call line K562 (calculated MW of 23.9kDa according to NP_001034750.1 and 24.8kDa according to NP_683708.1). Recommended concentration: 0.1-0.3ug/ml.	
Expiration Date	12 months from date of receipt.	