


Product Datasheet

STX6 Antibody (orb1246895)

Description	STX6 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-QALAERKNRQA	orb1246895 (0.3 ug/ml) staining of human...
Target	STX6	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 30kDa band observed in Human Lymph Node lysates (calculated MW of 29.2kDa according to NP_005810.1). In transfected HEK293 transiently expressing Human STX6 a band of approx. 30kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 0.3-1ug/ml.	
Clonality	Polyclonal	
Uniprot ID	O43752	
NCBI	NP_005810.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 30kDa band observed in Human Lymph Node lysates (calculated MW of 29.2kDa according to NP_005810.1). In transfected HEK293 transiently expressing Human STX6 a band of approx. 30kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 0.3-1ug/ml.	
Expiration Date	12 months from date of receipt.	