

Product Datasheet

STMN1 Antibody (orb1246902)

Description
STMN1 Antibody
Species/Host

Goat

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

EIA, ELISA, IHC, WB

Applications

orb1246902
(1 ug/ml)
staining of
Human T...

Immunogen

The immunogen for this antibody is: C-KNKESKDPADETED

Target

STMN1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage


Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:1000. Western Blot: Approx 17kDa band observed in Human Testis and Mouse Brain lysates (calculated MW of 17.3kDa according to NP_005554.1, NP_981944.1 and NP_981946.1). Recommended concentration: 1-3ug/ml. Immunohistochemistry: Paraffin embedded Human Thymus and Tonsil. Recommended concentration: 3.75ug/ml. ELISA: This antibody has been successfully used in ELISA as described in the following paper: Song et al, Aging Cell. 2010 Aug;9(4):607-15; PMID: 20560902.



orb1246902
(3.75 ug/ml)
staining of
para...



orb1246902
(3.75 ug/ml)
staining of
para...

Clonality

Polyclonal

Uniprot ID
P16949
NCBI
NP_005554.1, NP_981946.1, NP_981944.1
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:1000. Western Blot: Approx 17kDa band observed in Human Testis and Mouse Brain lysates (calculated MW of 17.3kDa according to NP_005554.1, NP_981944.1 and NP_981946.1). Recommended concentration: 1-3ug/ml. Immunohistochemistry: Paraffin embedded Human Thymus and Tonsil. Recommended concentration: 3.75ug/ml. ELISA: This antibody has been successfully used in ELISA as described in the following paper: Song et al, Aging Cell. 2010 Aug;9(4):607-15; PMID: 20560902.

Expiration Date

12 months from date of receipt.