



Product Datasheet

SMUG1 Antibody (orb1246938)

Description	SMUG1 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: PQAFLLSIHEPA-C	orb1246938 (0.1 ug/ml) staining of lysat...
Target	SMUG1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 38kDa band observed in lysates of cell line MOLT4 (calculated MW of 29.9kDa according to NP_055126.1). Recommended concentration: 0.1-0.3ug/ml. Primary incubation was 1 hour.	
Clonality	Polyclonal	
Uniprot ID	Q53HV7	
NCBI	NP_055126.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 38kDa band observed in lysates of cell line MOLT4 (calculated MW of 29.9kDa according to NP_055126.1). Recommended concentration: 0.1-0.3ug/ml. Primary incubation was 1 hour.	
Expiration Date	12 months from date of receipt.	



orb1246938
(5 ug/ml)
staining of
paraffi...