


Product Datasheet

TRIM40 Antibody (orb1247006)

Description	TRIM40 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-GVSELLLQPPQKL	orb1247006 (0.5 ug/ml) staining of Human...
Target	TRIM40	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 26kDa band observed in Human Placenta lysates (calculated MW of 26.1kDa according to NP_619645.1). Recommended concentration: 0.5-2ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q6P9F5	
NCBI	NP_619645.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 26kDa band observed in Human Placenta lysates (calculated MW of 26.1kDa according to NP_619645.1). Recommended concentration: 0.5-2ug/ml.	
Expiration Date	12 months from date of receipt.	