




Product Datasheet

RGS18 Antibody (orb1247021)

Description	RGS18 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	orb1247021 (1 ug/ml) staining of Human P...
Immunogen	The immunogen for this antibody is: CNEFQDVQSDVAIWL	
Target	RGS18	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247021 (5 ug/ml) staining of paraffi...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 28kDa band observed in Human Peripheral Blood Mononucleocytes (lymphocytes) lysates (calculated MW of 27.6kDa according to NP_570138.1). Recommended concentration: 1-3ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Prostate and Small Intestine. Recommended concentration: 5ug/ml.	
Clonality	Polyclonal	orb1247021 (5 ug/ml) staining of paraffi...
Uniprot ID	Q9NS28	
NCBI	NP_570138.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 28kDa band observed in Human Peripheral Blood Mononucleocytes (lymphocytes) lysates (calculated MW of 27.6kDa according to NP_570138.1). Recommended concentration: 1-3ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Prostate and Small Intestine. Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	