

Product Datasheet

PUM2 Antibody (orb1247050)

Description

PUM2 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

The immunogen for this antibody is: C-DQNRDLKQGDDDDSK

Target

PUM2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: In transfected HEK293 transiently expressing full-length Human PUM2 (myc and DYKDDDDK tagged), a band of approx. 140kDa was observed. No bands were observed in mock-transfected HEK293 and the same band was observed using anti-DYKDDDDK tag antibody. Recommended concentration 0.3-0.5ug/ml. Immunohistochemistry: Paraffin embedded Human Testis. Recommended concentration: 10ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q8TB72](#)

NCBI


[NP_056132.1](#)

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: In transfected HEK293 transiently expressing full-length Human PUM2 (myc and DYKDDDDK tagged), a band of approx. 140kDa was observed. No bands were observed in mock-transfected HEK293 and the same band was observed using anti-DYKDDDDK tag antibody. Recommended concentration 0.3-0.5ug/ml. Immunohistochemistry: Paraffin embedded Human Testis. Recommended concentration: 10ug/ml.

Expiration Date

12 months from date of receipt.



HEK293 lysate
(10 ug protein
in RIPA buf...



orb1247050
(10 ug/ml)
staining of
paraff...