




Product Datasheet

PRDM11 Antibody (orb1247077)

Description	PRDM11 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: LSEGRVRSGLCGG	orb1247077 (2 ug/ml) staining of A549 ly...
Target	PRDM11	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247077 (5 ug/ml) staining of paraffi...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 50kDa band observed in Human Lung Carcinoma A549 cell line lysates (calculated MW of 58kDa according to NP_064614.2). Recommended concentration: 2-3ug/ml. Immunohistochemistry: Paraffin embedded Human Brain (Cortex). Recommended concentration: 5ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q9NQV5	
NCBI	NP_064614.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 50kDa band observed in Human Lung Carcinoma A549 cell line lysates (calculated MW of 58kDa according to NP_064614.2). Recommended concentration: 2-3ug/ml. Immunohistochemistry: Paraffin embedded Human Brain (Cortex). Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	