




Product Datasheet

PIR Antibody (orb1247097)

Description	PIR Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	orb1247097 (0.3 ug/ml) staining of Human...
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: GSSKKVTLVSLSREC	
Target	PIR	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	orb1247097 (5 ug/ml) staining of paraffi...
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 33kDa band observed in Human Heart lysate (calculated MW of 32.1kDa according to NP_003653.1). Recommended concentration: 0.3 - 1.0 ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Liver and Spleen. Recommended concentration: 5ug/ml.	
Clonality	Polyclonal	orb1247097 (5 ug/ml) staining of paraffi...
Uniprot ID	O00625	
NCBI	NP_003653.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 33kDa band observed in Human Heart lysate (calculated MW of 32.1kDa according to NP_003653.1). Recommended concentration: 0.3 - 1.0 ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Liver and Spleen. Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	