

Product Datasheet

PSMF1 Antibody (orb1247105)

Description	PSMF1 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-DHLPPPGYDDMYL	orb1247105 staining (1 ug/ml) of Human K...
Target	PSMF1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 35kDa band observed in Human Kidney lysates (calculated MW of 29.8kDa according to NP_006805.2). In transfected HEK293 transiently expressing PSMF1 a band of approx. 35kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 1-3ug/ml.	HEK293 overexpressing PSMF1 and probed w...
Clonality	Polyclonal	
Uniprot ID	Q92530	
NCBI	NP_006805.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 35kDa band observed in Human Kidney lysates (calculated MW of 29.8kDa according to NP_006805.2). In transfected HEK293 transiently expressing PSMF1 a band of approx. 35kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 1-3ug/ml.	
Expiration Date	12 months from date of receipt.	