



Product Datasheet

PDLIM4 Antibody (orb1247123)

Description	PDLIM4 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: C-DVVAVYPNAKVELV	orb1247123 staining (0.2 ug/ml) of A431 ...
Target	PDLIM4	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247123 (3.75 ug/ml) staining of para...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 37kDa band observed in A431 lysates (predicted size of approx. 35kDa according to NP_003678). Recommended for use at 0.1-0.5ug/ml. Immunohistochemistry: Paraffin embedded Human Breast. Recommended concentration: 3.75ug/ml.	
Clonality	Polyclonal	
Uniprot ID	P50479	
NCBI	NP_003678.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 37kDa band observed in A431 lysates (predicted size of approx. 35kDa according to NP_003678). Recommended for use at 0.1-0.5ug/ml. Immunohistochemistry: Paraffin embedded Human Breast. Recommended concentration: 3.75ug/ml.	
Expiration Date	12 months from date of receipt.	