


Product Datasheet

PARP2 Antibody (orb1247144)


Description

PARP2 Antibody


Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The immunogen for this antibody is: C-LDLFEVEKDGEKE
Target	PARP2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 70kDa band observed in Human Spleen lysates (calculated MW of 60.8kDa according to NP_005475.1. Recommended concentration: 1-3ug/ml. Immunohistochemistry: Paraffin embedded Human Skin and Adrenal Gland. Recommended concentration: 5ug/ml.
Clonality	Polyclonal
Uniprot ID	Q9UGN5
NCBI	NP_005475.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 70kDa band observed in Human Spleen lysates (calculated MW of 60.8kDa according to NP_005475.1. Recommended concentration: 1-3ug/ml. Immunohistochemistry: Paraffin embedded Human Skin and Adrenal Gland. Recommended concentration: 5ug/ml.
Expiration Date	12 months from date of receipt.



orb1247144
(1 ug/ml)
staining of
Human S...



orb1247144
(5 ug/ml)
staining of
paraffi...



orb1247144
(5 ug/ml)
staining of
paraffi...