

Product Datasheet

PDZD2 Antibody (orb1247147)

Description

PDZD2 Antibody

Species/Host

Goat

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

The immunogen for this antibody is:
PITQDNAVHLHLPLC

Target

PDZD2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments gave no signal but low background in Hela and A431 extracts at up to 1ug/ml. Please note that we would not expect our routine western blot protocol to reliably detect proteins as large as the predicted size of this protein (approx 304kDa according to NP_055837) We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Clonality

Polyclonal

Uniprot ID

[O15018](#)

NCBI

[NP_835260.2](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments gave no signal but low background in Hela and A431 extracts at up to 1ug/ml. Please note that we would not expect our routine western blot protocol to reliably detect proteins as large as the predicted size of this protein (approx 304kDa according to NP_055837) We would appreciate any feedback from people in the field - have any results been reported with