

Product Datasheet

OSBPL3 Antibody (orb1247165)

Description

OSBPL3 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, IF

Applications
Immunogen

The immunogen for this antibody is: MSDEKNLGVSQKL-C

Target

OSBPL3

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of A431 cells. Recommended concentration: 10ug/ml.

Clonality



Polyclonal

Uniprot ID
[Q9H4L5](#)
NCBI
[NP_663162.1](#), [NP_663160.1](#), [NP_663161.1](#), [NP_056365.1](#)
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:32000. Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of A431 cells. Recommended concentration: 10ug/ml.

Expiration Date

12 months from date of receipt.


 orb1247165
Immunofluorescence
analysis o...

 orb1247165 Flow
cytometric analysis
of p...