




Product Datasheet

SNX1 Antibody (orb1247188)

Description	SNX1 Antibody	 250kDa 150kDa 100kDa
Species/Host	Goat	 50kDa 37kDa
Reactivity	Human	 25kDa 20kDa 15kDa
Conjugation	Unconjugated	orb1247188 (1 ug/ml) staining of HeLa ly...
Tested Applications	ELISA, ICC, IHC, WB	
Immunogen	The immunogen for this antibody is: C-TNSSKPQTYEELE	
Target	SNX1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	 250kDa 150kDa 100kDa 75kDa 50kDa 37kDa 25kDa 20kDa 15kDa
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247188 (3.8 ug/ml) staining of paraf...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 75kDa band observed in Human lysates of cell line HeLa (calculated MW of 59.1kDa according to NP_003090.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 1-3ug/ml. Immunohistochemistry: In paraffin embedded Human Uterus shows vesicular staining in the cytoplasm of endometrial epithelium. Recommended concentration, 3-5ug/ml. Immunocytochemistry: Vesicular staining of SNX1 (green) in formaldehyde-fixed SH5Y5Y. Recommended concentration, 3-5ug/ml.	 orb1247188 (5 ug/ml) staining of formald...
Clonality	Polyclonal	
Uniprot ID	Q13596	
NCBI	NP_690039.1 , NP_003090.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 75kDa band observed in Human lysates of cell line HeLa (calculated MW of 59.1kDa according to NP_003090.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 1-3ug/ml. Immunohistochemistry: In paraffin embedded Human Uterus shows vesicular staining in the cytoplasm of endometrial epithelium. Recommended concentration, 3-5ug/ml. Immunocytochemistry: Vesicular staining of SNX1	