

Product Datasheet

SLC9A3R2 Antibody (orb1247200)

Description

SLC9A3R2 Antibody



Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb1247200
(1 ug/ml)
staining of
Daudi l...

Immunogen

The immunogen for this antibody is: C-DLPGSDKDTEESGLH

Target

SLC9A3R2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 35kDa band observed in lysates of cell lines HEK293 and Daudi (calculated MW of 36.2kDa according to NP_004776.3). Recommended concentration: 1-3ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q15599](#)

NCBI

[NP_004776.3](#), [NP_001239005.1](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 35kDa band observed in lysates of cell lines HEK293 and Daudi (calculated MW of 36.2kDa according to NP_004776.3). Recommended concentration: 1-3ug/ml.

Expiration Date

12 months from date of receipt.