

## Product Datasheet

### MYRIP Antibody (orb1247231)

## Description

MYRIP Antibody

250kDa  
150kDa  
100kDa

### Species/Host

Goat

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, IHC, WB

### Applications

orb1247231  
(0.3 ug/ml)  
staining of  
Human...

### Immunogen

The immunogen for this antibody is: C-KDLMEPALESAVMY

### Target

MYRIP

### Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

### Form/Appearance

Liquid

### Concentration

500 ug/mL

### Storage


Aliquot and store at -20°C. Minimize freezing and thawing.

### Note


For research use only

### Application notes

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 90kDa band observed in Human Brain (Amygdala) lysates (calculated MW of 95.7kDa according to NP\_056275.2). Recommended concentration: 0.3-1ug/ml. Immunohistochemistry: Paraffin embedded Human Breast, Colon and Brain (Cortex). Recommended concentration: 3.75ug/ml.



orb1247231  
(3.75 ug/ml)  
staining of  
para...



### Clonality

Polyclonal

### Uniprot ID

[Q8NFW9](#)

### NCBI

[NP\\_056275.2](#)

orb1247231  
(3.75 ug/ml)  
staining of  
para...

### Dilution Range

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 90kDa band observed in Human Brain (Amygdala) lysates (calculated MW of 95.7kDa according to NP\_056275.2). Recommended concentration: 0.3-1ug/ml. Immunohistochemistry: Paraffin embedded Human Breast, Colon and Brain (Cortex). Recommended concentration: 3.75ug/ml.

### Expiration Date

12 months from date of receipt.