


Product Datasheet

MTR Antibody (orb1247244)

Description	MTR Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	The immunogen for this antibody is: C-VEKWLGPILGYDTD
Target	MTR
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000.Immunohistochemistry:In paraffin embedded Human Adrenal Gland shows cytoplasm staining in the cortex. Recommended concentration: 3-5ug/ml.
Clonality	Polyclonal
Uniprot ID	Q99707
NCBI	NP_000245.2 , NP_001278868.1 , NP_001278869.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000.Immunohistochemistry:In paraffin embedded Human Adrenal Gland shows cytoplasm staining in the cortex. Recommended concentration: 3-5ug/ml.
Expiration Date	12 months from date of receipt.



orb1247244
(3.8 ug/ml)
staining of
paraf...