



Product Datasheet

MTMR9 Antibody (orb1247245)

Description	MTMR9 Antibody	 
Species/Host	Goat	
Reactivity	Rat	orb1247245 (5 ug/ml) staining of paraffi...
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC	
Immunogen	The immunogen for this antibody is: C-AELETEDGMQESP	
Target	MTMR9	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: No signal obtained yet but low background observed in Human Brain and Kidney lysates at up to 1ug/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Immunohistochemistry: In paraffin embedded Rat Spleen shows staining of cytoplasm in select cells. Recommended concentration, 5-10ug/ml. This product has also been successfully used on Rat Hypothalamus.	
Clonality	Polyclonal	
Uniprot ID	Q96QG7	
NCBI	NP_056273.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: No signal obtained yet but low background observed in Human Brain and Kidney lysates at up to 1ug/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Immunohistochemistry: In paraffin embedded Rat Spleen shows staining of cytoplasm in select cells. Recommended concentration, 5-10ug/ml. This product has also been successfully used on Rat Hypothalamus.	
Expiration Date	12 months from date of receipt.	