

Product Datasheet

MORF4L1 Antibody (orb1247253)

Description

MORF4L1 Antibody

Species/Host

Goat

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

ELISA, WB
The immunogen for this antibody is: APKQDPKPKFQEG-C

Target

MORF4L1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 37kDa band observed in nuclear lysates of cell lines HepG2 and NIH3T3 and approx 40kDa band observed in Mouse, Rat and Pig Liver lysates (calculated MW of 37.2kDa according to Human NP_006782.1, Mouse NP_077751.1, Rat NP_001011999.1 and Pig XP_001927776.1 and 41.5kDa according to Human NP_996670.1, Mouse NP_001034236.1 and Pig XP_001927789.1). Recommended concentration: 0.1-0.3ug/ml.

Clonality

Polyclonal

Uniprot ID


[Q9UBU8](#)

NCBI


[NP_006782.1](#), [NP_996670.1](#)

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 37kDa band observed in nuclear lysates of cell lines HepG2 and NIH3T3 and approx 40kDa band observed in Mouse, Rat and Pig Liver lysates (calculated MW of 37.2kDa according to Human NP_006782.1, Mouse NP_077751.1, Rat NP_001011999.1 and Pig XP_001927776.1 and 41.5kDa according to Human NP_996670.1, Mouse NP_001034236.1 and Pig XP_001927789.1). Recommended concentration: 0.1-0.3ug/ml.



orb1247253
(0.3 ug/ml)
staining of
NIH3T...



orb1247253
(0.3 ug/ml)
staining of
Mouse...



orb1247253
(2.5 ug/ml)
staining of
paraf...