

Product Datasheet

MID2 Antibody (orb1247266)

Description

MID2 Antibody

Species/Host

Goat

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

ELISA, IHC, WB
The immunogen for this antibody is: QESPYVSGMKTCH

Target

MID2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 80-85kDa band observed in Mouse and Rat Heart, and Rat Kidney lysates (calculated MW of 77.8kDa according to Mouse NP_035975.1 and 83.1kDa according to Rat NP_001178818.1). Recommended concentration: 0.1-0.3ug/ml. Immunohistochemistry: In paraffin embedded Human Kidney shows strong staining in cytoplasm and brush border of some epithelial cells of proximal convoluted tubules. Recommended concentration, 5-10ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q9UJV3](#)

NCBI


[NP_036348.2](#), [NP_438112.2](#)

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 80-85kDa band observed in Mouse and Rat Heart, and Rat Kidney lysates (calculated MW of 77.8kDa according to Mouse NP_035975.1 and 83.1kDa according to Rat NP_001178818.1). Recommended concentration: 0.1-0.3ug/ml. Immunohistochemistry: In paraffin embedded Human Kidney shows strong staining in cytoplasm and brush border of some epithelial cells of proximal convoluted tubules. Recommended concentration, 5-10ug/ml.

Expiration Date

12 months from date of receipt.



orb1247266
(0.1 ug/ml)
staining of
Mouse...



orb1247266
(5 ug/ml)
staining of
paraffi...