

Product Datasheet

LC3A/B antibody (PE) (orb124727)

Description

Rabbit polyclonal antibody against LC3A/B

Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit, Rat
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human LC3A/B
Target	MAP1LC3B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Application notes	orb124727 is the PE conjugated form of orb6295
Isotype	IgG
Clonality	Polyclonal
MW	13kDa
Uniprot ID	Q9GZQ8
Expiration Date	12 months from date of receipt.