

Product Datasheet

MAN2A1 Antibody (orb1247289)

Description

MAN2A1 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

The immunogen for this antibody is: C-KISSDIKSQNR

Target

MAN2A1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx. 150kDa band observed in Human Placenta lysates (calculated MW of 131kDa according to NP_002363.2). In transfected HEK293 transiently expressing MAN2A1 a band of approx. 98kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 1-3ug/ml. An anonymous customer found positive results in WB on Human dermal fibroblasts.

Clonality

Polyclonal

Uniprot ID

[Q16706](#)

NCBI

[NP_002363.2](#)

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx. 150kDa band observed in Human Placenta lysates (calculated MW of 131kDa according to NP_002363.2). In transfected HEK293 transiently expressing MAN2A1 a band of approx. 98kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 1-3ug/ml. An anonymous customer found positive results in WB on Human dermal fibroblasts.

Expiration Date

12 months from date of receipt.




orb1247289
(0.3 ug/ml)
staining of
Human...



HEK293
overexpressing
Man2A1 and
probed ...