



Product Datasheet

MAML1 Antibody (orb1247292)

Description		MAML1 Antibody	
Species/Host	Goat		
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	ELISA, IHC, WB		
Immunogen	The immunogen for this antibody is: VLPTCPMAEFALPR		orb1247292 (0.1 ug/ml) staining of Human...
Target	MAML1		
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.		
Form/Appearance	Liquid		
Concentration	500 ug/mL		
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.		orb1247292 (3.8 ug/ml) staining of paraf...
Note	For research use only		
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000.Western Blot:Approx 110kDa band observed in Human Skeletal Muscle lysates (calculated MW of 108kDa according to NP_055572.1). Recommended concentration: 0.1-0.3ug/ml.Immunohistochemistry:In paraffin embedded Human Placenta shows nuclear staining in trophoblasts. Recommended concentration, 2.5-3.8ug/ml.		
Clonality	Polyclonal		
Uniprot ID	Q92585		
NCBI	NP_055572.1		
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000.Western Blot:Approx 110kDa band observed in Human Skeletal Muscle lysates (calculated MW of 108kDa according to NP_055572.1). Recommended concentration: 0.1-0.3ug/ml.Immunohistochemistry:In paraffin embedded Human Placenta shows nuclear staining in trophoblasts. Recommended concentration, 2.5-3.8ug/ml.		
Expiration Date	12 months from date of receipt.		