

Product Datasheet

Tspo Antibody (orb1247339)

Description

Tspo Antibody

Species/Host

Goat

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IF, WB

Applications
Immunogen

The immunogen for this antibody is: C-RDNSGRRGGSRLPE

Target

Tspo

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes


Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx. 17+18kDa doublet band observed in cell lysates of NIH3T3 and in Rat Adrenal Gland lysates (calculated MW of 18.8kDa according to Mouse NP_033905.3 and 18.9kDa according to Rat NP_036647.1). Recommended concentration: 1-2ug/ml. Primary incubation 1 hour at room temperature. Immunofluorescence: Strong expression of the protein seen in the cytoplasm of NIH3T3 and HeLa cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of NIH3T3 cells. Recommended concentration: 10ug/ml.

Clonality


Polyclonal

Uniprot ID
P50637
NCBI
NP_033905.3
Dilution Range


Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx. 17+18kDa doublet band observed in cell lysates of NIH3T3 and in Rat Adrenal Gland lysates (calculated MW of 18.8kDa according to Mouse NP_033905.3 and 18.9kDa according to Rat NP_036647.1). Recommended concentration: 1-2ug/ml. Primary incubation 1 hour at room temperature. Immunofluorescence: Strong expression of the protein seen in the cytoplasm of NIH3T3 and HeLa cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of NIH3T3 cells. Recommended concentration: 10ug/ml.



orb1247339
(3.75 ug/ml)
staining of
para...



orb1247339
(3.75 ug/ml)
staining of
para...



orb1247339
(3.75 ug/ml)
staining of
para...