

Product Datasheet

DBNL Antibody (orb1247397)

Description		DBNL Antibody	
Species/Host	Goat		
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	ELISA, WB		
Immunogen	The immunogen for this antibody is: AANLSRNGPALQE-C		orb1247397 (0.03 ug/ml) staining of Huma...
Target	DBNL		
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.		
Form/Appearance	Liquid		
Concentration	500 ug/mL		
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.		
Note	For research use only		
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 55kDa band observed in lysates of Human Peripheral Blood Mononucleocytes (calculated MW of 48.3kDa according to NP_054782.2). Recommended concentration: 0.03-0.1ug/ml. Primary incubation was 1 hour.		
Clonality	Polyclonal		
Uniprot ID	Q9UJU6		
NCBI	NP_001116428.1, NP_054782.2, NP_001014436.1		
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 55kDa band observed in lysates of Human Peripheral Blood Mononucleocytes (calculated MW of 48.3kDa according to NP_054782.2). Recommended concentration: 0.03-0.1ug/ml. Primary incubation was 1 hour.		
Expiration Date	12 months from date of receipt.		