

Product Datasheet

UBE2K Antibody (orb1247398)

Description

UBE2K Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

The immunogen for this antibody is: C-SWDVETATELLSN

Target

UBE2K

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:2000. Western Blot: Approx 24kDa band observed in Jurkat cell lysates (calculated MW of 22.4 kDa according to NP_005330). Recommended for use at 1-3 ug/ml. Immunohistochemistry: In paraffin embedded Human Brain shows pixulate cytoplasm staining in neurons of the Cortex. Recommended concentration, 1-2ug/ml.

Clonality

Polyclonal

Uniprot ID

P61086

NCBI

NP_005330.1, NP_001104582.1, NP_001104583.1

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:2000. Western Blot: Approx 24kDa band observed in Jurkat cell lysates (calculated MW of 22.4 kDa according to NP_005330). Recommended for use at 1-3 ug/ml. Immunohistochemistry: In paraffin embedded Human Brain shows pixulate cytoplasm staining in neurons of the Cortex. Recommended concentration, 1-2ug/ml.

Expiration Date

12 months from date of receipt.




250kDa
150kDa
100kDa



50kDa
37kDa
25kDa
20kDa
15kDa

orb1247398
staining (1.5
ug/ml) of
Jurka...



orb1247398
(3.8 ug/ml)
staining of
paraf...