




Product Datasheet

NDC80 Antibody (orb1247402)

Description	NDC80 Antibody	 250kDa 150kDa 100kDa
Species/Host	Goat	 50kDa 37kDa 25kDa
Reactivity	Human	orb1247402 staining (0.5 ug/ml) of Hela
Conjugation	Unconjugated	...
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: C-YEKKATLIKSSEE	
Target	NDC80	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247402 (5 ug/ml) staining of paraffi...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 75-80kDa band observed in Hela lysate (predicted MW of 74kDa according to NP_006092). Recommended for use at 0.5-2ug/ml. Immunohistochemistry: Paraffin embedded Human Liver and Adrenal Gland. Recommended concentration: 5ug/ml.	
Clonality	Polyclonal	
Uniprot ID	O14777	orb1247402 (5 ug/ml) staining of paraffi...
NCBI	NP_006092.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 75-80kDa band observed in Hela lysate (predicted MW of 74kDa according to NP_006092). Recommended for use at 0.5-2ug/ml. Immunohistochemistry: Paraffin embedded Human Liver and Adrenal Gland. Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	