


Product Datasheet

GRIPAP1 Antibody (orb1247417)


Description

GRIPAP1 Antibody

Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	The immunogen for this antibody is: LREQTTGLAAELQ-C
Target	GRIPAP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments gave no signal but low background in Human Brain lysates at up to 1ug/ml. Immunohistochemistry: Paraffin embedded Human Tonsil and Brain (Cortex). Recommended concentration: 3.75ug/ml.
Clonality	Polyclonal
Uniprot ID	Q4V328
NCBI	NP_997555.1 , NP_064522.3
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments gave no signal but low background in Human Brain lysates at up to 1ug/ml. Immunohistochemistry: Paraffin embedded Human Tonsil and Brain (Cortex). Recommended concentration: 3.75ug/ml.
Expiration Date	12 months from date of receipt.



orb1247417
(3.75 ug/ml)
staining of
para...



orb1247417
(3.75 ug/ml)
staining of
para...