


Product Datasheet

GALR3 Antibody (orb1247464)

Description

GALR3 Antibody

Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The immunogen for this antibody is: C-QGPEPREGPVHGGE
Target	GALR3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 37kDa band observed in Human Skin and Human Testes lysates (calculated MW of 39.6kDa according to NP_003605.1). Recommended concentration: 1-3ug/ml. Immunohistochemistry: In paraffin embedded Human Testis shows cellular membrane staining of spermatocytes. Recommended concentration, 4-6ug/ml.
Clonality	Polyclonal
Uniprot ID	O60755
NCBI	NP_003605.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 37kDa band observed in Human Skin and Human Testes lysates (calculated MW of 39.6kDa according to NP_003605.1). Recommended concentration: 1-3ug/ml. Immunohistochemistry: In paraffin embedded Human Testis shows cellular membrane staining of spermatocytes. Recommended concentration, 4-6ug/ml.
Expiration Date	12 months from date of receipt.



orb1247464
(2 ug/ml)
staining of
Human S...



orb1247464
(3.8 ug/ml)
staining of
paraf...