



Product Datasheet

GAD2 Antibody (orb1247466)

Description	GAD2 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-TLEDNEERMSRLSK	orb1247466 (0.5 ug/ml) staining of Human...
Target	GAD2	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247466 (2.5 ug/ml) staining of paraf...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 65kDa band observed in Human Brain (Hippocampus) lysates (calculated MW of 65.4kDa according to NP_000809.1). Recommended concentration: 0.5-1.5ug/ml. Primary incubation was 1 hour.	
Clonality	Polyclonal	
Uniprot ID	Q05329	
NCBI	NP_000809.1 , NP_001127838.1	orb1247466 (2.5 ug/ml) staining of paraf...
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 65kDa band observed in Human Brain (Hippocampus) lysates (calculated MW of 65.4kDa according to NP_000809.1). Recommended concentration: 0.5-1.5ug/ml. Primary incubation was 1 hour.	
Expiration Date	12 months from date of receipt.	