



Product Datasheet

FBXL3 Antibody (orb1247522)

Description	FBXL3 Antibody	
Species/Host	Goat	 
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: KRGRDSDRNSSEE-C	
Target	FBXL3	orb1247522 (1 ug/ml) staining of Human A...
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 50-55kDa band observed in Human Adipose lysates (calculated MW of 48.7kDa according to NP_036290.1). Recommended concentration: 1-3ug/ml. Immunohistochemistry: Paraffin embedded Human Thyroid. Recommended concentration: 5ug/ml.	orb1247522 (5 ug/ml) staining of paraffi...
Clonality	Polyclonal	
Uniprot ID	Q9UKT7	
NCBI	NP_036290.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 50-55kDa band observed in Human Adipose lysates (calculated MW of 48.7kDa according to NP_036290.1). Recommended concentration: 1-3ug/ml. Immunohistochemistry: Paraffin embedded Human Thyroid. Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	