



Product Datasheet

FNBP1 Antibody (orb1247525)

Description	FNBP1 Antibody	 250kDa 150kDa
Species/Host	Goat	 75kDa 50kDa 37kDa
Reactivity	Human	 25kDa 20kDa
Conjugation	Unconjugated	orb1247525 staining (0.25 ug/ml) of Huma...
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: KQLESSKRRFERDC	
Target	FNBP1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247525 (3.75 ug/ml) staining of para...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 70kDa band observed in Human Testis lysates (calculated MW of 71.3kDa according to NP_055848.1). Recommended concentration: 0.25-2ug/ml. Immunohistochemistry: Paraffin embedded Human Tonsil and Small Intestine. Recommended concentration: 3.75ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q96RU3	orb1247525 (3.75 ug/ml) staining of para...
NCBI	NP_055848.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 70kDa band observed in Human Testis lysates (calculated MW of 71.3kDa according to NP_055848.1). Recommended concentration: 0.25-2ug/ml. Immunohistochemistry: Paraffin embedded Human Tonsil and Small Intestine. Recommended concentration: 3.75ug/ml.	
Expiration Date	12 months from date of receipt.	