


Product Datasheet

FBXL2 Antibody (orb1247526)

Description	FBXL2 Antibody	
Species/Host	Goat	
Reactivity	Human, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: VFSNNDEGLINKK-C	
Target	FBXL2	orb1247526 (0.3 ug/ml) staining of Human...
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 48kDa band observed in Human Brain (Frontal Cortex) lysates and approx 50kDa band observed in Rat Brain lysates (calculated MW of 47.1kDa according to Human NP_036289.2, and 46.9kDa according to Rat XP_343496.4). Recommended concentration: 0.3-1ug/ml. Primary incubation was 1 hour. Immunohistochemistry: In paraffin embedded Human Cerebral Cortex shows cytoplasm staining in neuronal bodies. Recommended concentration: 5-10ug/ml.	orb1247526 (10 ug/ml) staining of paraff...
Clonality	Polyclonal	
Uniprot ID	Q9UKC9	
NCBI	NP_036289.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 48kDa band observed in Human Brain (Frontal Cortex) lysates and approx 50kDa band observed in Rat Brain lysates (calculated MW of 47.1kDa according to Human NP_036289.2, and 46.9kDa according to Rat XP_343496.4). Recommended concentration: 0.3-1ug/ml. Primary incubation was 1 hour. Immunohistochemistry: In paraffin embedded Human Cerebral Cortex shows cytoplasm staining in neuronal bodies. Recommended concentration: 5-10ug/ml.	
Expiration Date	12 months from date of receipt.	