


Product Datasheet

XPO7 Antibody (orb1247539)

Description	XPO7 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	orb1247539 staining (0.5 ug/ml) of Hela ...
Immunogen	The immunogen for this antibody is: C-NSTYGVNSNDMMMS	
Target	XPO7	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247539 (5 ug/ml) staining of paraffi...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 110kDa band observed in Hela lysates (predicted size of approx. 124kDa according to NP_055839). Recommended concentration: 0.5-1ug/ml Immunohistochemistry: Paraffin embedded Human Testis. Recommended concentration: 5ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q9UIA9	
NCBI	NP_001093632.1 , NP_055839 , NP_001093631.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 110kDa band observed in Hela lysates (predicted size of approx. 124kDa according to NP_055839). Recommended concentration: 0.5-1ug/ml Immunohistochemistry: Paraffin embedded Human Testis. Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	