

Product Datasheet

FBXO43 Antibody (orb1247559)

Description

FBXO43 Antibody

Species/Host

Goat

Reactivity

Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb1247559
(0.3 ug/ml)
staining of
Mouse...

Immunogen

The immunogen for this antibody is: C-PGSAQSKRNLKRL

Target

FBXO43

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 80kDa band observed in Mouse Ovary lysates (calculated MW of 78.4kDa according to NP_001025031.2). Recommended concentration: 0.3-1ug/ml. Primary incubation was 1 hour. Preliminary testing was unsuccessful on Human, Rat and Pig Testis, and Human Ovary for this particular batch.

Clonality

Polyclonal

Uniprot ID

[Q4G163](#)

NCBI

[NP_001025031.2](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 80kDa band observed in Mouse Ovary lysates (calculated MW of 78.4kDa according to NP_001025031.2). Recommended concentration: 0.3-1ug/ml. Primary incubation was 1 hour. Preliminary testing was unsuccessful on Human, Rat and Pig Testis, and Human Ovary for this particular batch.

Expiration Date

12 months from date of receipt.

