

Product Datasheet

ELF1 Antibody (orb1247566)

Description

ELF1 Antibody

Species/Host

Goat

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

ELISA, IHC, WB
The immunogen for this antibody is: C-AMKQNELLEPNSF

Target

ELF1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 65kDa band observed in Human Lymph Node and Tonsil and in Mouse Thymus lysates (calculated MW of 64.9kDa according to Human NP_001138825.1 and 62.3kDa according to Mouse NP_001273341.1). Recommended concentration: 0.1-0.3ug/ml. Approx 70kDa band observed in lysates of cell line Jurkat (calculated MW of 67.5kDa according to NP_758961.1). Recommended concentration: 0.03-0.1ug/ml. Immunohistochemistry: Paraffin embedded Human Thymus. Recommended concentration: 2.5ug/ml.

Clonality

Polyclonal

Uniprot ID


[P32519](#)

NCBI


[NP_001138825.1](#), [NP_758961.1](#)

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 65kDa band observed in Human Lymph Node and Tonsil and in Mouse Thymus lysates (calculated MW of 64.9kDa according to Human NP_001138825.1 and 62.3kDa according to Mouse NP_001273341.1). Recommended concentration: 0.1-0.3ug/ml. Approx 70kDa band observed in lysates of cell line Jurkat (calculated MW of 67.5kDa according to NP_758961.1). Recommended concentration: 0.03-0.1ug/ml. Immunohistochemistry: Paraffin embedded Human Thymus. Recommended concentration: 2.5ug/ml.



orb1247566
(0.3 ug/ml)
staining of
Human...



orb1247566
(0.1 ug/ml)
staining of
Jurka...



orb1247566
(2.5 ug/ml)
staining of
paraf...