

Product Datasheet

DLC1 Antibody (orb1247600)

Description

DLC1 Antibody

Species/Host

Goat

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, IF, IHC

Applications

Immunogen

The immunogen for this antibody is: SVAIRKRSWEEHC

Target

DLC1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:64000. Immunohistochemistry: Paraffin embedded Mouse Adrenal Gland. Recommended concentration: 4ug/ml. Immunofluorescence: Expression of the protein seen in the cytoplasm and nucleus of U251 and U2OS cells. Recommended concentration: 10ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q96QB1](#)

NCBI



[NP_872584.2](#), [NP_079043.3](#)

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:64000. Immunohistochemistry: Paraffin embedded Mouse Adrenal Gland. Recommended concentration: 4ug/ml. Immunofluorescence: Expression of the protein seen in the cytoplasm and nucleus of U251 and U2OS cells. Recommended concentration: 10ug/ml.

Expiration Date


12 months from date of receipt.



orb1247600
Immunofluorescence
analysis o...



orb1247600
Immunofluorescence
analysis o...



orb1247600 (4
ug/ml) staining of
paraffi...