



Product Datasheet

Cyld Antibody (orb1247635)

Description	Cyld Antibody	
Species/Host	Goat	
Reactivity	Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: C-DSQQSKSKNPWYIDE	
Target	Cyld	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 125kDa band observed in Rat Brain lysates (calculated MW of 107kDa according to mouse NP_775545.1). Recommended concentration: 0.5-1.5ug/ml. Immunohistochemistry: Paraffin embedded Mouse Brain. Recommended concentration: 5ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q80TQ2	
NCBI	NP_775545.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 125kDa band observed in Rat Brain lysates (calculated MW of 107kDa according to mouse NP_775545.1). Recommended concentration: 0.5-1.5ug/ml. Immunohistochemistry: Paraffin embedded Mouse Brain. Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	