


Product Datasheet

CNTFR Antibody (orb1247664)

Description	CNTFR Antibody	
Species/Host	Goat	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: CRSNTYPKGFY	orb1247664 (0.3 ug/ml) staining of mouse...
Target	CNTFR	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 38-40kDa band observed in mouse brain lysates (calculated MW of 40.8kDa according to mouse NP_057882.1). Recommended concentration: 0.3-1ug/ml.	
Clonality	Polyclonal	
Uniprot ID	P26992	
NCBI	NP_001833.1 , NP_671693.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 38-40kDa band observed in mouse brain lysates (calculated MW of 40.8kDa according to mouse NP_057882.1). Recommended concentration: 0.3-1ug/ml.	
Expiration Date	12 months from date of receipt.	