


Product Datasheet

CILP Antibody (orb1247670)

Description		CILP Antibody	
Species/Host	Goat		
Reactivity	Rat		
Conjugation	Unconjugated		
Tested Applications	ELISA, WB		
Immunogen	The immunogen for this antibody is: C-NYRRTDHEDPRVK		orb1247670 (0.03 ug/ml) staining of rat ...
Target	CILP		
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.		
Form/Appearance	Liquid		
Concentration	500 ug/mL		
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.		
Note	For research use only		
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 150kDa band observed in rat spinal cord lysates (calculated MW of 133kDa according to human NP_003604.2 and of 132kDa according to Rat NP_001101631.1). Recommended concentration: 0.1-0.3ug/ml.		
Clonality	Polyclonal		
Uniprot ID	O75339		
NCBI	NP_003604.2		
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 150kDa band observed in rat spinal cord lysates (calculated MW of 133kDa according to human NP_003604.2 and of 132kDa according to Rat NP_001101631.1). Recommended concentration: 0.1-0.3ug/ml.		
Expiration Date	12 months from date of receipt.		