


Product Datasheet

BPNT1 Antibody (orb1247726)

Description	BPNT1 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-RNYDYYASRVESIK	orb1247726 (1 ug/ml) staining of Human L...
Target	BPNT1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Approx 37kDa band observed in Human Liver lysates (calculated MW of 33.4kDa according to NP_006076.4). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Spiegelberg et al, J Biol Chem. 1999 May 7;274(19):13619-28.; PMID: 10224133). Recommended concentration: 1-3ug/ml.	
Clonality	Polyclonal	
Uniprot ID	O95861	
NCBI	NP_006076.4	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Approx 37kDa band observed in Human Liver lysates (calculated MW of 33.4kDa according to NP_006076.4). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Spiegelberg et al, J Biol Chem. 1999 May 7;274(19):13619-28.; PMID: 10224133). Recommended concentration: 1-3ug/ml.	
Expiration Date	12 months from date of receipt.	