

Product Datasheet

RNF115 Antibody (orb1247737)

Description		RNF115 Antibody	
Species/Host	Goat		
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	ELISA, WB		orb1247737 (0.3 ug/ml) staining of HEPG2...
Immunogen	The immunogen for this antibody is: C-DQDNRRANERGHQTH		
Target	RNF115		
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.		
Form/Appearance	Liquid		
Concentration	500 ug/mL		
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.		
Note	For research use only		
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 35kDa band observed in Human lysates of cell line HEPG2 (calculated MW of 33.7kDa according to NP_055270.1). Recommended concentration: 0.3-1ug/ml.		
Clonality	Polyclonal		
Uniprot ID	Q9Y4L5		
NCBI	NP_055270.1		
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 35kDa band observed in Human lysates of cell line HEPG2 (calculated MW of 33.7kDa according to NP_055270.1). Recommended concentration: 0.3-1ug/ml.		
Expiration Date	12 months from date of receipt.		