

## Product Datasheet

### PPP2R5B Antibody (orb1247747)

## Description

PPP2R5B Antibody

### Species/Host

Goat

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

orb1247747  
(2 ug/ml)  
staining of  
Human C...

### Immunogen

The immunogen for this antibody is: C-QRLTPQVAASGGQS

### Target

PPP2R5B

### Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

### Form/Appearance

Liquid

### Concentration

500 ug/mL

### Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

### Note

For research use only

### Application notes

Peptide ELISA: antibody detection limit dilution 1:2000. Western Blot: Approx. 55kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 57.4kDa according to NP\_006235.1). Recommended concentration: 2-3ug/ml. Primary incubation was 1 hour. Preliminary testing was unsuccessful on Mouse and Rat Brain for this particular batch.

### Clonality

Polyclonal

### Uniprot ID

[Q15173](#)

### NCBI

[NP\\_006235.1](#)

### Dilution Range

Peptide ELISA: antibody detection limit dilution 1:2000. Western Blot: Approx. 55kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 57.4kDa according to NP\_006235.1). Recommended concentration: 2-3ug/ml. Primary incubation was 1 hour. Preliminary testing was unsuccessful on Mouse and Rat Brain for this particular batch.

### Expiration Date

12 months from date of receipt.

