


Product Datasheet

ARPC1B Antibody (orb1247764)

Description	ARPC1B Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	orb1247764 staining (2 ug/ml) of Human L...
Immunogen	The immunogen for this antibody is: AYHSFLVEPISCH	
Target	ARPC1B	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247764 (5 ug/ml) staining of paraffi...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 37kDa band observed in Human Liver lysates (calculated MW of 40.9kDa according to NP_005711). Recommended concentration: 2-4ug/ml. A minor band of unknown identity was also consistently observed at 28kDa. Immunohistochemistry: Paraffin embedded Human Thymus. Recommended concentration: 5ug/ml.	
Clonality	Polyclonal	
Uniprot ID	O15143	
NCBI	NP_005711.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 37kDa band observed in Human Liver lysates (calculated MW of 40.9kDa according to NP_005711). Recommended concentration: 2-4ug/ml. A minor band of unknown identity was also consistently observed at 28kDa. Immunohistochemistry: Paraffin embedded Human Thymus. Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	