


Product Datasheet

MKI67 antibody (orb12478)

Description	Rabbit polyclonal antibody to MKI67
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	A synthetic peptide of human Ki67(NP_001139438.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	359kDa
Uniprot ID	P46013
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...