




Product Datasheet

AKAP9 Antibody (orb1247818)

Description	AKAP9 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IF, WB	
Immunogen	The immunogen for this antibody is: C-SGSTTQFHAGMR	orb1247818 Immunofluorescence analysis o...
Target	AKAP9	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247818 Immunofluorescence analysis o...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: This product has been successfully used by a customer, showing a band at approx 450kDa in lysates of cell line HepG2 (calculated MW of 453kDa according to NP_005742.4). Recommended concentration: 1-3ug/ml. Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 and U2OS cells. Recommended concentration: 10ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q99996	
NCBI	NP_671714.1 , NP_005742.4	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: This product has been successfully used by a customer, showing a band at approx 450kDa in lysates of cell line HepG2 (calculated MW of 453kDa according to NP_005742.4). Recommended concentration: 1-3ug/ml. Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 and U2OS cells. Recommended concentration: 10ug/ml.	
Expiration Date	12 months from date of receipt.	