

Product Datasheet

ADRA2A Antibody (orb1247838)

Description

ADRA2A Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, FC, IF

Immunogen

The immunogen for this antibody is: C-TERRP NGLGPERS

Target

ADRA2A

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Immunofluorescence: Strong expression of the protein seen in the membranes of MCF7 cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of MCF7 cells. Recommended concentration: 10ug/ml.

Clonality

Polyclonal

Uniprot ID

[P08913](#)

NCBI


[NP_000672.3](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:32000. Immunofluorescence: Strong expression of the protein seen in the membranes of MCF7 cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of MCF7 cells. Recommended concentration: 10ug/ml.


Expiration Date

12 months from date of receipt.



orb1247838

Immunofluorescence analysis o...



orb1247838 Flow cytometric analysis of p...